

Sample Preparation for Sorting on the FACS Aria

| | | |
|------------------------------------|--|--|
| All Runs Require: | <ul style="list-style-type: none"> • All the samples must be in a single cell suspension and must have been filtered through a mesh of less than 70 microns. • Samples are preferred in sterilized PBS with 1% FBS/PBS. • To avoid the clumping, DNase can be added (50ug/ml) in the final cell suspension. • Include 7AAD/PI to stain your dead cells (if your viability <75%). • For optimum sorting and save your precious samples, a diluted aliquot (300-400ul total) of your cells-to-be sorted should be placed in a separate tube in order to set up the instrument. | |
| Controls: | <u>Negative unstained controls</u> (0.3ml minimum volume) | <u>Compensation controls for each fluor</u> (0.3ml minimum volume) |
| Acquisition Samples Volume: | Cells for sorting should be in 0.5ml minimum volume if you have less than 1×10^6 cells. | Data acquisition/Record runs typically only acquire 10,000~30,000 gated events. |
| Sort Cell Concentration: | Samples should be at $10 \sim 20 \times 10^6$ cells/ml for high speed sorts and $5 \sim 10 \times 10^6$ /ml for low and medium speed sorts (5,000 to 35,000 events/s) | The Aria can sort at up to 20~100 million cells/hour (a flow rate of 0.2-0.8ul/s and up to 2ul/s with less resolution). |
| Collection Tubes: | Sorted cells can be accepted in 1 ml, 5 ml, or 15 ml tubes. 15ml tubes must be BD or Corning. | The tubes should contain medium with antibiotics and 50% of FBS if healthy and viable sorted cells are preferred. |
| Attention: | For the first time such as particular cell sorting or preparing cells, it is best to plan a trial run of your experiment. | Be aware that if the samples are not prepared properly, then problems such as clumping and bad staining may occur; your sort may not be completed within the reserved time; and your flow profiles may not meet your expectations. |
| Questions: | Please direct all questions to Lisa Xu at (604) 875-2000 ext 5987 <Lixin86@interchange.ubc.ca > | |